




# SUPER CONFORT SET

XT Auto *Neo*  
HYBERNITE

XT Auto Neoスーパーコンフォートセット

**乾燥した冬やアレルギーの季節も快適なCPAP治療を！  
CPAP治療コンプライアンス向上のために！**


XT Auto専用の加温加湿器とハイバーナイト加温回路により、常に最適な加温加湿状態を保ち、回路内の結露を防ぎます。



使用者が“主観的に最も快適と感じる”<sup>\*1</sup> 加温加湿状態で治療を行うことができ、CPAP治療のコンプライアンス向上がはかれます。

従来の加温加湿器のみ使用したCPAP治療と比べ、加温回路により加湿状態がコントロールされた場合、深睡眠、レム睡眠ならびに総睡眠時間が増加したという報告<sup>\*2</sup>や鼻腔内抵抗が低減した結果CPAPタイトレーション圧が平均10%低減したという報告<sup>\*3</sup>、鼻腔内ステロイドと鼻の症状改善に対して同等の効果があったという報告<sup>\*4</sup>もあります。

CPAP使用中の諸症状に対するハイバーナイトシステムの使用に関して、  
下記のようなアンケート調査結果が出ております\*5。



\*1 "Climate control: humidification with heated tube", A V Benjafield et al, ResMed Science Center, ResMed Ltd, Sydney, Australia  
 \*2 "Impact of a controlled heated breathing tube humidifier on sleep quality during CPAP therapy in a cool sleeping environment", G Niluis et al, European Respiratory Journal 2008; 31: 830-836  
 \*3 "Effect of humidification on titration pressure in obstructive sleep apnea", J S Massengill et al, Sleep 2009; 32(suppl): A217  
 \*4 "Heated breathing tube vs nasal steroids for compliance and quality of life", E Powell et al, Sleep 2010; 33(suppl): A469&A483  
 \*5 Survey on Hybernite usage for clinical symptoms during CPAP therapy. PlastiFlex, Jan/Feb 2011



総販売元

株式会社 MAGnet

〒160-0022 東京都新宿区新宿 6-27-49 新宿風越ビル

TEL: 03-3200-7033 FAX: 03-3200-7034

<http://www.magnet-japan.com/>